

Myocardial Protection *At the heart of the Protection*



Dedicated solutions for all of the different myocardial protection methods. A wide collection of crystalloid, cold and warm blood cardioplegia systems available for delivering your myocardial protection.

1. Cold Crystalloid Cardioplegia

Developed in the mid -1960s, cold crystalloid cardioplegia¹ represented a key milestone in open heart surgery, providing a simple but effective way of protecting the heart: hypothermia. In many institutions around the world, cold crystalloid cardioplegia remains in use as a method of myocardial protection.

LivaNova offers several sets for general use, and customizes sets upon request.

1. H. Geissler, U. Mehlhorn "Cold Crystalloid Cardiopneuma", multimedia Cardio Thoracic Surgery volume 2006, Issue 0109.

Order Guide		
Item Number	Product Designation	Units per box
AB1610	D500 - Basic Clear Cardioplegia Delivery Set	10
AB1611	D505 - Optimized Clear Cardioplegia Delivery Set	2

2. Blood Cardioplegia

Developed as an evolution of cold crystalloid cardioplegia, blood cardioplegia (cold or warm) provides the following benefits:

- Improved oxygen carrying capacity and delivery
- Enhanced myocardial oxygen consumption
- Preserved high-energy phosphate stores
- Buffering changes in pH
- · Provides appropriate osmotic environment for myocardial cells and lessens the myocardial oedema

	Oxygen Demand	Reduction ²
Normothermic Arrest (37°C)	1ml/100g/min	90%
Hypothermic Arrest (22°C) Set	0.30 ml/100g/min	97%
Hypothermic Arrest (10°C)Set	0.14 ml/100g/min	~97%

2.1. Warm Blood Cardioplegia with Bubble Trap

The D725 Warm Blood Cardioplegia delivery set has been developed in response to the Modern Protocols targeting myocardial protection by lowering the oxygen consumption due to chemical stop (-90% in normothermia)².



The dedicated BCD Vanguard[™] bubble trap improves safety for the patient ameliorating the air handling.

The BCD Vanguard bubble trap includes a 105 µm screen filters. Its special geometry filters and directs bubbles up and out of the device through a hydrophobic membrane and a one-way valve. The result is enhanced safety and the convenience of one step priming. A pressure relief valve (> +400 mmHg) prevents any damage generated by occlusion of the table line.

2. GD Buckberg, JR Brazier and RL Nelson et al., Studies of the effects of hypothermia in regional blood flow and metabolism during cardiopulmonary bypass, J Thorac Cardiovasc Surg 73 (1977), pp. 87-94.

BCD Bubble Trap

Technical Features

Connections	3/16" (inlet & outlet)
Static Priming Volume - BCD bubble trap only	23 ml
Static Priming Volume - D725 circuit	130 ml
Air Handling	- integrated bubble trap - 105 μm filter screen
Max Blood Flow	500 ml/min

Order Guide		
Item Number	Product Designation	Units per box
05095	D725 Warm Blood Cardioplegia Delivery Set	15

2.2. Cold Blood Cardioplegia: CSC and Vanguard

When the combination of hypothermia security and optimum oxygen delivery are needed, cold blood cardioplegia is the response.

LivaNova offers a wide range of different stainless steel heat exchangers to fully meet your needs.

> CSC 14



The Visible Solution

Total visibility design is a common concept in the LivaNova ECC offering. The clear fluid pathways increase the possibility for overall visible monitoring by perfusionists, thereby enhancing safety during the extracorporeal perfusion.

Like LivaNova oxygenators, the LivaNova Venous and Cardiotomy Reservoirs as well as the LivaNova Arterial Filters, the LivaNova Cardioplegia Heat Exchangers are designed with this philosophy in mind: the visible solution.

Great performance

The CSC 14 Heat Exchanger is designed to have an extremely high performance factor for single pass perfusion both for blood cardioplegia and crystalloid deliveries.

Efficient and quick preparation

The priming and debubbling of CSC 14 is very quick and easy, enhanced by the integrated priming dial and complete visibility of the fluid pathway.

Integrated Bubble Trap

An air separation screen creates a highly efficient bubble trap, minimizing the risk of air bubble delivery.

Unique Reversible Holder

The holder allows right or left pump set-up and is affixed from the back side of CSC 14, maintaining easy access and full visibility.

Flexibility and Versatility

CSC 14 is available separate or integrated in standard cardioplegia sets.

1.0 0.8 0.6 0.4 0.2 0 500 1000 1500 2000 2500 FLOW RATE (ML/MIN.)

HEAT EXCHANGE PERFORMANCE

Technical Features

Connections	- 1/4″ Blood inlet - 1/4″ Blood Outlet - 1/2″ Hansen Quick
Static Priming Volume - Device only	30 ml
Pressure drop @ 300 ml/min	5 mm Hg
Air Handling	- integrated Bubble trap - 150 μm filter screen
Housing Material	Polycarbonate
Max Blood Flow	2500 ml/min

Order Guide		
Item Number	Product Designation	Units per box
095808	CSC 14 Cardioplegia Set 1:1	10
C20173	CSC 14 Cardioplegia Set 4:1	5
C20174	CSC 14 Cardioplegia Set 4:1 with Shunt	5
C20175	CSC 14 Cardioplegia Set 4:1 (pediatric)	5
C20176	CSC 14 Cardioplegia Set 1:1 (pediatric)	5
050137	HSC Holder	5
P3740	CSC 14 sterile stand alone	5

> BCD Vanguard™ Efficiency, Safety and Convenience



Efficiency

The BCD Vanguard[™] Heat Exchanger incorporates a chevron-grooved, closely packed, pleated stainless steel sheet that enhances mixing without hemolytic turbulence.

Thanks to the limited sheet and blood film thickness, the heat transfer is highly efficient, making hypothermic arrest and cardioplegic rewarming extremely rapid.

Safety

Safety is strongly emphasized throughout the design of the BCD Vanguard[™] heat exchanger. Its pressure relief valve prevents rupture if a table line is inadvertently clamped, helping protect both the device and the patient.

The device has a bubble trap that includes a filter screen to deflect bubbles up to the membrane, which then vents the air out of the system.



Convenience

The hydrophobic membrane of the BCD Vanguard[™] makes priming more convenient without any extra steps. By clamping the outlet line and turning on the pump, the air is automatically purged out of the device via the membrane. This unique feature reduces unnecessary delays.





Cardiopulmonary Bypass patients may also benefit from tip-to-tip **Phosphorylcholine Coating**



P.H.I.S.I.O. coated Vanguard[™] blood cardioplegia sets project patients decreases blood activation to patients heart during cardioplegia perfusion.



Technical Features

Connections	- 1/4″ Blood inlet - 3/16″ Blood Outlet - 1/2″Hansen Quick Connectors (water ports)
Static Priming Volume - Device only	39 ml
Pressure drop @ 250 ml/min	<10 mm Hg
Air Handling	- integrated Bubble trap - 105 µm filter screen
Housing Material	Polycarbonate
Max Blood Flow	500 ml/min

Order Guide		
Item Number	Product Designation	Units per box
05471	D921 BCD Vanguard 1:1 Cardioplegia Set	5
05472	D924 Vanguard 4:1 Cardioplegia Set	5
05473	D924S Vanguard 4:1 Cardioplegia Set with shunt	5
05474	D924P Vanguard 4:1 Cardioplegia Set pediatric	5
05475	D928S Vanguard 8:1 Cardioplegia Set with shunt	5
C21243	Fully coated BCD Vanguard™ 4:1 Cardioplegia Set	5
050137	HSC holder	

www.livanova.com



Advanted by

The Sorin Group Italia Quality System complies with: EN ISO 13485:2012 C C 0123 According to Annex II (Ful Quality System) of MDD 93/42/EFC as amended by directive 2007/47/EFC



© 2017 LivaNova all rights reserved.

Manufactured by: Sorin Group Italia S.r.l. a wholly-owned subsidiary of LivaNova PLC Via Statale 12 Nord, 86 41037 Mirandola (MO) Italia Tel.: +39.0535.29811 info.cardiacsurgery@livanova.com

