SORIN

3T AEROSOL COLLECTION SET

en - ENGLISH - INSTRUCTIONS FOR USE.....	3	cs - ČEŠTINA - NÁVOD K POUŽITÍ.....	51
it - ITALIANO - ISTRUZIONI PER L'USO	7	pl - POLSKI - INSTRUKCJA OBSŁUGI	55
fr - FRANÇAIS - MODE D'EMPLOI.....	11	sk - SLOVENČINA - NÁVOD NA POUŽITIE	59
de - DEUTSCH - GEBRAUCHSANWEISUNG	15	sl - SLOVENSKO - NAVODILA ZA UPORABO	63
es - ESPAÑOL - INSTRUCCIONES DE USO	19	hu - MAGYAR - KEZELÉSI ÚTMUTATÓ	67
pt - PORTUGUÊS - INSTRUÇÕES DE UTILIZAÇÃO.....	23	lt - LIETUVIŠKAI - NAUDOJIMO INSTRUKCIJOS	71
el - ΕΛΛΗΝΙΚΑ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	27	et - EESTI - KASUTUSJUHEND.....	75
nl - NEDERLANDS - GEBRUIKSAANWIJZINGEN	31	lv - LATVISKI - LIETOŠANAS INSTRUKCIJAS	79
sv - SVENSKA - BRUKSANVISNING	35	ro - ROMÂNĂ - INSTRUCȚIUNI DE UTILIZARE.....	83
da - DANSK - BRUGSANVISNING	39	bg - БЪЛГАРСКИ - ИНСТРУКЦИИ ЗА УПОТРЕБА	87
fi - SUOMI - KÄYTTÖOHJEET.....	43	tr - TÜRKÇE - KULLANMA TALİMATI.....	91
no - NORSK - BRUKSANVISNING.....	47		

3T AEROSOL COLLECTION SET

REF: 050900100

This document supplements the current IFU for the Heater-Cooler System 3T and should be used in conjunction with that IFU.

INTRODUCTION

The 3T Aerosol Collection Set collects aerosol created by the Heater-Cooler 3T water tanks during use. This 3T Aerosol Collection Set replaces the use of the overflow bottle referenced in the Heater-Cooler System 3T Operating Instructions. **Please read all information carefully:** Failure to properly follow these instructions may allow for aerosol to be emitted from the Heater-Cooler 3T in its surrounding environment.

INDICATIONS FOR USE

The Stöckert Heater-Cooler System 3T is used with a Stöckert S3 heart-lung machine and/or any other heart-lung machine featuring a separate temperature control for extracorporeal perfusion of durations up to 6 hours.

DEVICE

The 3T Aerosol Collection Set consists of a suction canister, a suction canister lid and two pieces of tubing; the longer piece of tubing includes a connector (see Figure 1).


Figure 1) Overview of the 3T Aerosol Collection Set: 1 = suction canister, 2 = pre-connected right angle connector, 3 = Heater-Cooler connection line, 4 = vacuum source line

WARNINGS AND PRECAUTIONS

- The use of a central hospital vacuum source is recommended for use with this disposable. The vacuum source must be capable of 20 l/min flow or greater.
- If a central hospital vacuum source is not available the use of a portable vacuum system capable of 20 l/min flow per ISO 10079-1 for Medical Suction Equipment is required. The portable vacuum source must be specified for medical applications and fulfill the requirements of IEC 60601-1 and IEC 60601-1-2.
- Do not place the portable vacuum system or any other object on top of the Heater-Cooler 3T.
- Follow the portable vacuum system supplier's mounting and operating instructions for the recommended use of the pump in an operating room environment.
- To avoid disturbances in the laminar flow area of the operating room, make sure that a portable vacuum system is placed in such a way that the exhaust flow is not directed towards the operating field, but towards the exhaust vent system.
- Power sources for the vacuum system selected must have connection to a back-up power generator to avoid loss of power during use
- If the Patient #2 circuit is used and water tank level adjustment is required, only add pre-mixed water/H₂O₂ solution as per the separate Heater-Cooler 3T operating instructions if the water level drops to the level where the red LED for the patient tank is activated. In such a case, add water until the first green LED is visible.
- The 3T Aerosol Collection Set is validated for 3 calendar days of use.
- Check the date of installation, which is written on the suction canister lid as per point 13 of this "INSTRUCTIONS FOR USE", for use-by-date of the 3T Aerosol Collection Set.
- The suction canister is not intended as a measuring device.
- The suction canister contents are potentially hazardous due to the possible collection of aerosolized bacteria. Use appropriate personal protection equipment (PPE) and handle accordingly.
- Store in dark place: long term exposure to light may compromise product performance and result in breakage during use.

INSTRUCTIONS FOR USE

WARNING indicates serious adverse reactions and potential safety hazards for the practitioner and/or the patient that may occur in the proper use or misuse of the device as well as the limitations of use and the measures to be adopted in such cases.

CAUTION indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

1. Before starting a procedure, with the heater-cooler and the vacuum source switched off, verify that the tank vacuum gauge indicator located on the Heater-Cooler is positioned at 0 (see Figure 2). If necessary, set the indicator needle to 0 using the adjustment screw near the bottom of the gauge.

Figure 2) Tank Vacuum gauge with zero point and adjustment screw in the bottom part of the gauge

2. Ensure that a suction canister mounting bracket is attached to the right side or rear panel of the Heater-Cooler 3T. The mounting bracket is attached by inserting the two top tabs into the vent slots on the Heater-Cooler 3T and then allowing the bracket to rest in a horizontal position. The mounting bracket should be positioned so that the base of the suction canister is above the base of the Heater-Cooler 3T to avoid displacement of the canister from the holder during use. Please refer to the diagrams below.

Figure 3) Suction canister Mounting bracket

Figure 4) Mounting bracket installed on right side panel.

Figure 5) Mounting bracket installed on rear panel.


Figure 6) Tilt mounting bracket up to insert or remove.

3. Remove the components from the packaging.
4. Place the lid on the canister and press it down firmly around entire perimeter (see Figure 7).

Figure 7) Canister lid assembled on the canister


Figure 8) Overview of the canister lid components: P = Heater-cooler connection port, V = vacuum port, T = tandem port (shown here with tandem port cap applied), S = pour spout (shown here with pour spout cap applied), L = label for recording the installation date. The hydrophobic filter is located inside the device below the vacuum port.

5. Apply the cap firmly over the matching pour spout (S) (see Figures 7 and 8).
6. Apply the tandem port cap tightly over the tandem port (T) (see Figures 7 and 8).
7. Place the suction canister with the mounted lid into the suction canister holder (see Figures 9 and 10).


Figure 9) Canister with lid inserted into the suction canister holder on side panel

Figure 10) Canister with lid inserted into the suction canister holder on rear panel


Figure 11) Overview of the 3T Aerosol Collection Set: 1 = suction canister, 2 = pre-connected right angle connector, 3 = Heater-Cooler connection line, 4 = vacuum source line

8. Attach the short 46 cm (1.5 ft) Heater-Cooler connection tube 3 (see Figure 11) to the Heater-Cooler Overflow Outlet and then connect it to the port (P) of the canister lid (see Figures 8, 9 and 10).
Caution: Do not connect the short 46 cm (1.5 ft) Heater-Cooler connection tube to the center vacuum source port (V) (see Figure 8) as this will potentially cause the hydrophobic filter to be exposed to water and will limit the vacuum flow to the canister.
9. Attach the 366 cm (12 ft) 1/4" ID vacuum source line 4 with pre-connected right angle connector 2 (see Figure 11) to vacuum port (V) (see Figures 8, 9 and 10), making sure that the tubing is fully connected and no tubing is kinked.
10. If additional length is required to access the vacuum source, attach the 366 cm (12 ft) 1/4" ID vacuum extension line with connector to the 366 cm (12 ft) 1/4" ID vacuum source line 4 (see Figure 11).
Note: If required, the vacuum extension line must be ordered separately; REF 050900111.
11. Attach the other end of the 366 cm (12 ft) 1/4" ID vacuum source line 4 (see Figure 11) to a regulated vacuum source (see Figure 12).

Caution: The vacuum source line must be connected to a vacuum regulator to control and adjust flow to the 3T Aerosol Collection Set. Do not exceed vacuum level of 640 mm Hg (25" Hg).

Note: It is recommended to use an overflow safety trap connected to the regulated vacuum source to avoid liquid introduction into the medical–surgical vacuum system.


Figure 12) 3T Aerosol Collection Set connected to the 3T and to the regulated vacuum source with the overflow safety trap

12. Check all caps and connections for proper sealing.
13. Write the date of installation on the label (L) located on the lid (see Figure 8)
Caution: The 3T Aerosol Collection Set must be replaced within 3 days of installation. It must also be replaced after a Heater-Cooler 3T water change or disinfection procedure.

Figure 13) Tank Vacuum gauge with indicator in the acceptable area of operation and adjustment screw in the bottom part of the gauge

14. With the Heater-Cooler 3T switched off, adjust the vacuum source so that the Tank Vacuum gauge on the Heater-Cooler 3T is adjusted to approximately 50 Pa (see Figure 13).
15. Switch on the Heater-Cooler 3T and during use adjust the vacuum source to ensure that the Tank Vacuum gauge on the Heater-Cooler 3T remains between 21 Pa and 100 Pa (see Figure 13).

Warning: Operating the Heater-Cooler 3T with the vacuum level indicated on the Tank Vacuum gauge in the red solid zone limits the disposable's ability to collect aerosol from the Heater-Cooler 3T water tanks.

Refer to section "Troubleshooting" for possible remedies.

Note: It is advised to add the check of vacuum level indicated on the Tank Vacuum gauge to your procedure set-up check list.

16. Do not disconnect the vacuum source until the Heater-Cooler 3T is switched off. When the Heater-Cooler 3T is switched off, disconnect the vacuum source line from the vacuum regulator.

Warning: Operation of the Heater-Cooler 3T without vacuum source applied will stop collection of aerosol from the water tanks.

17. At the end of the procedure in order to prevent an excessive overflow of the tank volume capacity, empty the external circuits according to the following procedure:

Warning: Do not perform the external circuit emptying steps until the patient has been discharged from the Operating Room environment to eliminate the potential of exposure to aerosol.

- Ensure that the vacuum source is applied and the 3T Aerosol Collection Set is connected
- Close the Cardioplegia circuit valve and empty the Cardioplegia circuit as per the separate Heater-Cooler 3T operating instructions (water will move back to the Cardioplegia Tanks)
- Close the Patient #1 circuit valve and empty the Patient #1 circuit as per the separate Heater-Cooler 3T operating instructions (water will move back to the Patient Tank and potentially start overflowing from the tank to the canister)
- If the Patient #2 circuit is used and the Patient #1 circuit emptying has already filled the suction canister with a volume greater than 0.5 liter, proceed as described in the following steps:

- Empty the suction canister per the following steps:

Warning: As the canister contents are potentially hazardous due to collection of aerosolized bacteria, use appropriate Personal Protection Equipment (PPE) and handle the contents accordingly.

Warning: Do not expose the suction canister hydrophobic lid filter to the water inside the canister as this condition will seal the filter stopping vacuum flow and limiting the disposable's ability to collect aerosol from the Heater-Cooler 3T water tanks.

Warning: Do not lift canister by lid as the weight of the contents may cause the lid to separate from the canister

- Switch off the heater-cooler 3T.
- Disconnect the short 46 cm (1.5 ft) Heater-Cooler connection tube 3 (see Figure 11) from the port (P) of the suction canister lid (see Figures 8, 9 and 10)
- Disconnect the 366 cm (12 ft) 1/4" ID vacuum source line 4 with pre-connected right angle connector 2 (see Figure 11), from the vacuum port (V) of the suction canister lid (see Figures 8, 9 and 10)
- Open the pour spout cap (S) (see Figures 7 and 8)
- Empty the suction canister content according to hospital policy
- Connect again the short 46 cm (1.5 ft) Heater-cooler connection tube 3 (see Figure 11) to the port (P) of the suction canister lid (see Figures 8, 9 and 10)
- Connect again the 366 cm (12 ft) 1/4" ID vacuum source line 4 with pre-connected right angle connector 2 (see Figure 11), to the vacuum port (V) of the suction canister lid (see Figures 8, 9 and 10)
- Apply the cap firmly over matching pour spout (S) (see Figures 7 and 8).
- Switch on the heater-cooler 3T to allow empty of the Patient #2 circuit.

- Or drain water from the Patient tank through the tank drain valve of the Patient circuit until the orange LED blinks; once the orange LED is blinking, close the tank drain valve
- Close circuit valve and empty the Patient #2 circuit as per the separate Heater-Cooler 3T operating instructions (water will move back to the Patient Tank).

18. If during external circuit emptying the suction canister is filled with water from the Heater-Cooler Overflow Outlet, empty the suction canister following the procedure as described in item 17 above.

Caution: If the suction canister is completely filled with water and has exposed the hydrophobic filter located below the vacuum source port to water, replace the 3T Aerosol Collection set to avoid loss of vacuum performance caused by obstruction of the hydrophobic filter contacting the water. Exposure of the hydrophobic filter to water will cause the filter to block flow of air and will reduce and stop vacuum flow to the canister, reducing the disposable's ability to collect aerosol from the Heater-Cooler 3T water tanks.

DISPOSAL

Dispose of the canister and tubing within 3 calendar days of installation. Check the installation date which is written on the label on the suction canister lid as discussed in step 13 of "INSTRUCTIONS FOR USE".

Caution: When performing the Heater-Cooler 3T disinfection procedure leave the used 3T Aerosol Collection Set attached without vacuum source applied.

To avoid disinfectant vapor getting to the vacuum source, disconnect the 366 cm (12 ft) 1/4" ID vacuum source line from the suction canister prior to performing the disinfection procedure.

Replace the 3T Aerosol Collection Set after the disinfection procedure has been completed.

1. Switch off the Heater-Cooler 3T.
2. Disconnect the short 46 cm (1.5 ft) Heater-cooler connection tube 3 (see Figure 11) from the Heater-cooler Overflow Outlet (see Figures 9 and 10).
3. Disconnect the 366 cm (12 ft) 1/4" ID vacuum source line 4 (see Figure 11) from the vacuum source (see Figure 12).
4. Remove the suction canister from the suction canister holder.
5. Transport the 3T Aerosol Collection Set to the disposal area.
Warning: do not lift canister by lid as the weight of the contents may cause the lid to separate from the canister
6. Dispose the 3T Aerosol Collection Set according to hospital policy.

TROUBLESHOOTING

- If the measured pressure inside the Heater-Cooler 3T tank is not greater than 21 Pa (see Figure 13), adjust the vacuum regulator of the vacuum source until the Tank Vacuum gauge indicator reaches the acceptable area of operation (see Figure 13).
- If there is a loss of vacuum during the procedure check that the vacuum source is operational, the canister is properly sealed, that all connections are tight and tubing is not kinked.
- If the suction canister is completely filled with water, exchange the 3T Aerosol Collection set to avoid loss of vacuum performance caused by the obstruction of the hydrophobic filter contacting the water. Exposure of the hydrophobic filter to water will cause the filter to block flow of air and will reduce and stop vacuum flow to the canister, reducing the disposable's ability to collect aerosol from the 3T water tanks.
- If there is a loss of vacuum during the procedure the air flow to the vacuum source may be obstructed. Check the vacuum source and the overflow safety trap for possible obstruction.
- If no obstruction is found and a loss of vacuum is still observed, open the vacuum gauge service port located on the back of the machine (see Figure 14). Opening this service port by removing the luer lock cap will allow vacuum source air flow to remove any potential obstruction within the pressure monitoring line internal to the 3T. Re-attach the luer lock cap and confirm the obstruction has been removed.


Figure 14) Vacuum Gauge Service Port

- If loss of vacuum is still observed apply positive air pressure to the service port using a sterile syringe to deliver approximately 30cc of air into the service port. Re-attach the luer lock cap and confirm the obstruction has been removed.

PART NUMBERS

3T Aerosol Collection Set	050900100
¼" ID vacuum extension line with connector, 366 cm (12 ft)	050900111

The 3T Aerosol Collection Set is supplied non-sterile for 3 calendar days of use. Discard in accordance with hospital procedures after use.

LABELING

	Attention, see instruction for use
	Follow instructions for use (white symbol on blue background)
	Batch Code
	Purchase order number
Rx Only	Only applies in the U.S.A.: Sale (and prescription) is restricted to physicians
EA	Quantity
	Manufacturer
	Use By

	Dispose the 3T Aerosol Collection Set within 3 calendar days from the installation date written on the label
	Keep away from heat
	Keep dry

RETURN OF USED PRODUCT

For Customers within the United States

If for any reason this product must be returned to Sorin Group USA, Inc., a returned goods authorization (RGA) number is required from Sorin Group USA, Inc., prior to shipping.

If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton, or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of the shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from the Sorin Group USA, Inc. Returned Goods Coordinator, Quality Assurance Department (800-650-2623).

CAUTION: It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

The shipping address for returned goods is:

Sorin Group USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599

For Customers outside the United States

If for any reason this product must be returned, please contact your sales representative for specific instructions.

If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton, or an equivalent carton, to prevent damage during shipment.

CAUTION: It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

3T Aerosol Collection Set • Operating instructions

Copyright © 2016
SORIN GROUP DEUTSCHLAND GmbH Lindberghstrasse 25
D-80939 Munich, Germany

Tel.: +49/(0)89/32301-0
Fax: +49/(0)89/32301-555

All rights reserved, especially the right to reproduction and distribution as well as translation. No part of this document may be reproduced – by photocopy, microfilm or any other process – nor may any part of it be stored, edited, duplicated or distributed by electronic means without the written permission of SORIN GROUP DEUTSCHLAND GmbH.

The name "Stöckert®" is a Registered Trademark.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Manufactured for:
SORIN GROUP DEUTSCHLAND GmbH Lindberghstrasse 25
D-80939 Munich, Germany
Tel.:+49/(0)89/32301-0
Fax:+49/(0)89/32301-555

Distributed in the U.S.A. by: Sorin Group USA, Inc.
14401 W. 65th Way, Arvada, CO 80004,
phone:800.221.7943 phone:303.425.5508
fax:303.467.6584

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Copyright© 2016
SORIN GROUP USA, Inc